

GREAT LAKES ORTHOPEDIC LABS, INC.

2362 Genesee Street, Cheektowaga, NY 14225

Telephone. No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.

219 Bryant Street, Buffalo, NY 14222

Telephone No. (716) 878-7307

**INSTRUCTIONS FOR WEARING THE
ABDUCTION PILLOW**

1. Child should always have clothes on between pillow and skin.
2. Lie child supine.
3. Spread legs and insert pillow carefully between legs.
4. Attach each Velcro strap – around the outside of the thigh and attach to the front of the pillow.
5. **WASHING INSTRUCTIONS** – Clean orthosis with wash cloth and everyday bar soap two or three times per week.
6. Wear the brace the number of hours per day as prescribed by the doctor. Also follow physician's orders with respect to activities allowed.
8. **Caution:** When not wearing the device, keep it away from your pets since they are attracted to the human scent and will chew and destroy it.
9. Call immediately for an appointment if there are any skin problems, discomfort or other fitting problems. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.

Warranty: The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.