

GREAT LAKES ORTHOPEDIC LABS, INC.

2362 Genesee Street, Cheektowaga, NY 14225

Telephone. No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.

219 Bryant Street, Buffalo, NY 14222

Telephone No. (716) 878-7307

**INSTRUCTIONS FOR WEARING THE
CAM WALKER**

1. Patient should wear a heavy cotton sock free of wrinkles.
2. While sitting place your foot in the heel of the cam walker.
3. Secure the front panel inside (outside) of the brace.
4. Secure the straps from the bottom to the top making sure the straps are snug and comfortable.
5. **WASHING INSTRUCTIONS** - Wash the inside of the brace everyday with mild soap and water. Rinse and towel dry.
6. Wear the brace the number of hours per day as prescribed by the doctor. Also follow the doctor's orders with respect to activities allowed.
7. **Caution:** When not wearing the device, keep it away from your pets since they are attracted to the human scent and will chew and destroy it.
7. Call immediately for an appointment if there are any skin problems, discomfort or other fitting problems. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.

Warranty: The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.