

GREAT LAKES ORTHOPEDIC LABS, INC.
2362 Genesee Street, Cheektowaga, NY 14225
Telephone. No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.
219 Bryant Street, Buffalo, NY 14222
Telephone No. (716) 878-7307

INSTRUCTIONS FOR WEARING KNEE IMMOBILIZERS

1. Identify the top and bottom of the knee immobilizer.
2. Patient should be on his or her back.
3. Place knee immobilizer behind the leg and align with the knee.
4. Overlap the flaps of the knee immobilizer. The straps should be across from each other.
5. Check the alignment of the brace.
6. First tighten the middle straps, then the bottom straps, and finally the top strap. Tighten all straps once again as needed.
7. **WASHING INSTRUCTIONS** – Wash by hand with a gentle detergent. Drip dry.
8. Wear the brace the number of hours per day as prescribed by the doctor. Also follow the doctor's orders with respect to activities allowed.
9. Call immediately for an appointment if there are any skin problems, discomfort or other fitting problems. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.
10. **Caution:** When not wearing the device, keep it away from your pets since they are attracted to the human scent and will chew and destroy it.

Warranty: The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.