

**GREAT LAKES ORTHOPEDIC LABS, INC.**

2362 Genesee Street, Cheektowaga, NY 14225

Telephone. No. (716) 893-4116

**GREAT LAKES ORTHOPEDIC LABS, INC.**

219 Bryant Street, Buffalo, NY 14222

Telephone No. (716) 878-7307

**DIRECTIONS FOR WEARING THE  
PAVLIK HARNESS**

1. Open up the chest strap and place your baby on the chest strap at nipple level.
2. Fasten the chest strap with the velcro closure on the front of the chest. Make sure this strap is not too tight. You should be able to slide two fingers under the strap once it is tightened.
3. Attach the shoulder straps to the loops on the front chest strap.
4. Place the feet in the booties so that one strap is behind the knee and one strap is in front of the knee.
5. Make sure the legs are at 90° to the body and spread apart as directed by your doctor.
6. Keep the Pavlik Harness away from pets. Pets will destroy the Pavlik Harness if they chew on it.
7. As your baby grows, the shoulder, chest, and leg straps may need to be lengthened to accommodate growth.
8. If you are unsure of the correct position of your baby's legs or of the fit of the Pavlik Harness, please contact your doctor or call our office to schedule an appointment.
9. If you are having fitting problems, please call our office. If you need immediate assistance after regular office hours, please call (716) 893-4116 and leave a message on the answering machine. Your call will be returned.

**Warranty:** The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.