

GREAT LAKES ORTHOPEDIC LABS, INC.
2362 Genesee Street, Cheektowaga, NY 14225
Telephone No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.
219 Bryant Street, Buffalo, NY 14222
Telephone No. (716) 878-7307

INSTRUCTIONS FOR WEARING A STUMP SHRINKER

1. Roll shrinker in half.
2. Grip the shrinker and pull over the stump until bottom of shrinker encapsulates the bottom of the stump making sure there are no wrinkles.
3. Unroll rest of shrinker.
4. Shrinker should feel snug.
5. **WASHING INSTRUCTIONS:** Wash in cold water with mild soap. Rinse thoroughly under running water, removing all soapy residue. Drip dry.
6. Wear the shrinker as prescribed by the doctor. Also follow doctor's orders with respect to follow-up appointments.
7. Call immediately for an appointment if there are any skin problems, discomfort or other fitting problems. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.

Warranty: The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.