

GREAT LAKES ORTHOPEDIC LABS, INC.

2362 Genesee Street, Cheektowaga, NY 14225

Telephone No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.

219 Bryant Street, Buffalo, NY 14222

Telephone No. (716) 878-7307

DIRECTIONS FOR WEARING THUMB SPLINTS

- 1) Place the thumb in the neoprene thumb section of the splint.
- 2) Slide the thumb section towards the palm as far as possible so that the tip of the thumb sticks out of the thumb section.
- 3) Fasten the velcro strap on the back of the hand snugly so the thumb section stays in position close to the palm.
- 4) Follow the manufacturer's washing instructions for hygiene.
- 5) If you have fitting problems call our office for an appointment. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.
- 6) Caution: When not wearing the device, keep it away from your pets since they are attracted to the human scent and will chew and destroy it.

Warranty: The components of this device, since it is a prefabricated device, carry an individual manufacturer's warranty. Great Lakes will honor warranties under applicable law. Most devices are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.