

GREAT LAKES ORTHOPEDIC LABS, INC.

2362 Genesee Street, Cheektowaga, NY 14225

Telephone. No. (716) 893-4116

GREAT LAKES ORTHOPEDIC LABS, INC.

219 Bryant Street, Buffalo, NY 14222

Telephone No. (716) 878-7307

**DIRECTIONS FOR WEARING
DYNAMIC (ULTRAFLEX) ELBOW SPLINTS**

1. Place upper arm section on front of arm first.
2. Work the thumb into the thumb hole section of the wrist section.
3. Fasten the elbow pad.
4. Fasten the wrist/hand strap.
5. Wash the brace with soap and water using soap such as ivory or other soaps safe for skin use.
6. Caution: When not wearing the device, keep it away from your pets since they are attracted to the human scent and will chew and destroy it.
7. Monitor for skin pressure problems. All pressure marks should disappear within twenty minutes, once the brace is removed. If pressure marks last longer than twenty minutes please call for an appointment so that adjustments can be made to relieve the pressure areas.
8. If you are having any problems call 878-7307 for our office in Women & Children's Hospital or 893-4116 for our Cheektowaga office. If you need immediate assistance after regular office hours, please call 893-4116 and leave a message on the answering machine. Your call will be returned.

Warranty: The components of this device are warranted for 90 days, after which a service charge, based on an hourly rate and materials, can be applied. The warranty does not apply to anatomical changes of the patient, misuse and/or loss of the device or any alterations made by anyone other than Great Lakes.